

Sample Information


Name : 乙酰基四肽-3
Sequence : Ac-KGHK-NH₂ (Acetate Salt)
Pump A : 0.1% Trifluoroacetic in 100% water
Pump B : 0.1% Trifluoroacetic in 100% acetonitrile
Total Flow : 1ml/min
Wavelength : 214nm
Analytical column type : Kromasil 100-5-C18 (4.6*250mm*5 µm)

Dissolution method : 100% H₂O

Inj. Volume : 20 µl

Time	Module	Action	Value
0.01	Pumps	B.Conc	0.1
30.00	Pumps	B.Conc	30
33.00	Pumps	B.Conc	100
38.00	Pumps	B.Conc	100
40.00	Pumps	B.Conc	0.1
50.00	Controller	Stop	

Chromatogram



Peak Table

Detector A 214nm

Peak#	Ret. Time	Area	Height	Area%
1	4.956	31162	4688	0.184
2	7.690	239499	21881	1.412
3	8.324	16589728	1436830	97.804
4	8.865	101866	4091	0.601
Total		16962255	1467491	100.000

邮箱 : marketing@motifbiotech.com

网址 : <https://www.motifbiotech.com>

地址 : 苏州市高新区浒墅关镇浒青路 111 号 B04