

CERTIFICATE OF ANALYSIS

Product Name	乙酰基四肽-5
Sequen	Ac-β-{Ala}-HSH
Dissolution condition	100%H ₂ O
Length	4AA
Modification	acetate salt
Molecular Weight (MW)	492.50
Storage	-20°C

Test Items	Specifications	Results
MW by MS	493.25(Peptide)	Conforms
Purity by HPLC	>98%	98.70%
Peptide content	≤ 7.0%	1.77%
Moisture content	N/A	N/A
Acetic acid content	≤ 20.0%	13.77%
Appearance	White to off-white powd	Conforms
Quantity	3g	3.0g

Certified by:
Quality Assurance Department



Note: this product is intended for research use only; not for diagnostic or human use

邮箱: marketing@motifbiotech.com 网址: <https://www.motifbiotech.com>

地址: 苏州市高新区浒墅关镇浒青路 111 号 B04