

CERTIFICATE OF ANALYSIS

Product Name	十肽-12
Sequence	YRSRKYSSWY (Acetate Salt)
Dissolution condition	100%H2O
Length	10AA
Modification	acetate salt
Molecular Weight (MW)	1395.52
Storage	-20°C

Test Items	Specifications	Results
MW by MS	1395.30	Conforms
Purity by HPLC	>95%	98.926%
Peptide content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	1g	50.0mg*20管

Certified by:

Quality Assurance Department



Note: this product is intended for research use only; not for diagnostic or human use.

邮箱：marketing@motifbiotech.com

网址：https://www.motifbiotech.com

地址：苏州市高新区浒墅关镇浒青路 111 号 B04