

CERTIFICATE OF ANALYSIS

Product Name	二肽-2
Sequence	VW
Dissolution condition	100%H ₂ O
Length	2AA
Modification	sodium salt
Molecular Weight (MW)	303.36
Storage	-20°C

Test Items	Specifications	Results
MW by MS	303.34(Peptide)	Conforms
Purity by HPLC	>98%	99.54%
Peptide content	≥ 80%	88.71%
Moisture content	≤ 7.0%	3.35%
Acetic acid content	N/A	N/A
Appearance	White to off-white powd	Conforms
Quantity	2g	2.0g

Certified by:
Quality Assurance Department



Note: this product is intended for research use only; not for diagnostic or human use.

邮箱: marketing@motifbiotech.com 网址: <https://www.motifbiotech.com>

地址: 苏州市高新区浒墅关镇浒青路 111 号 B04