

CERTIFICATE OF ANALYSIS

Product Name	棕榈酰四肽-3
Sequence	Pal-GQPR (HCl Salt)
Dissolution condition	5%HCOOH+15%ACN+80%H2O
Length	4AA
Modification	HCl Salt
Molecular Weight (MW)	694.91
Storage	-20°C

Test Items	Specifications	Results
MW by MS	694.40	Conforms
Purity by HPLC	>95%	98.004%
Peptide content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	1g	50.0mg*20管

Certified by:

Quality Assu

**Note: this product is intended for research use only; not for diagnostic or human use.**

邮箱：marketing@motifbiotech.com

网址：https://www.motifbiotech.com

地址：苏州市高新区浒墅关镇浒青路 111 号 B04