

CERTIFICATE OF ANALYSIS

Product Name	五肽-3
Sequence	GPRPA-NH2 (Acetate Salt)
Dissolution condition	100%H2O
Length	5AA
Modification	Acetate salt
Molecular Weight (MW)	495.56
Storage	-20°C

Test Items	Specifications	Results
MW by MS	495.60	Conforms
Purity by HPLC	>95%	97.649%
Peptide Content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	40mg	5.0mg*8管

Certified by:

Quality Assurance Department

**Note: this product is intended for research use only; not for diagnostic or human use.**网址: <https://www.motifbiotech.com>邮箱: marketing@motifbiotech.com

地址: 苏州市高新区浒墅关镇浒青111号 B04