




Sample Information

Dissolution method	:15%ACN+85%H ₂ O	Interface	:ESI	Prerod Bias	:+4.5kv
Date Acquired	:2023/12/12 15:06:24	Nebulizing Gas Flow	:1.50L/min	Detector	:-0.2kv
Injection Volume	:0.2ul	CDL Temp	:250°C	T.Flow	:0.2ml/min
Block Temp	:200	CDL Volt	:0v	B.conc	:50%H ₂ O/50%MEOH

Name :五肽-3
 Sequence :GPRPA-NH₂ (Acetate Salt)
 Theoretical :495.56
 Observed :495.60

网址：<https://www.motifbiotech.com>

邮箱：marketing@motifbiotech.com

地址：苏州市高新区浒墅关镇浒青路111号B04