




Sample Information

Dissolution method :15%ACN+85%H<sub>2</sub>O  
Date Acquired :2024/04/23 16:35:48  
Injection Volume :0.2ul  
Block Temp :200

Interface :ESI  
Nebulizing Gas Flow :1.50L/min  
CDL Temp :250°C  
CDL Volt :0v

Prerod Bias :+4.5kv  
Detector :-0.2kv  
T.Flow :0.2ml/min  
B.conc :50%H<sub>2</sub>O/50%MEOH

Name :六肽-9  
Sequence :GPQGPQ  
Theoretical :582.62  
Observed :582.40

网址：<https://www.motifbiotech.com>

邮箱：[marketing@motifbiotech.com](mailto:marketing@motifbiotech.com)

地址：苏州市高新区浒墅关镇浒青路111号B04