


Sample Information

Name : 蓝铜胜肽
Sequence : GHK-Cu (Acetate Salt)
Pump A : 0.1% Trifluoroacetic in 100% water
Pump B : 0.1% Trifluoroacetic in 100% acetonitrile
Total Flow : 1ml/min
Wavelength : 214nm
Analytical column type : YMC-Triart C18(4.6*250mm*5um)
Dissolution method : 100%H2O
Inj. Volume : 45ul

Time	Module	Action	Value
0.01	Pumps	B.Conc	0.1
20.00	Pumps	B.Conc	0.1
33.00	Pumps	B.Conc	100
38.00	Pumps	B.Conc	100
40.00	Pumps	B.Conc	0.1
50.00	Controller	Stop	

Chromatogram



Peak Table

Detector A 214nm

Peak#	Ret. Time	Area	Height	Area%
1	3.609	9131331	725505	98.185
2	4.063	160577	3191	1.727
3	8.560	8193	811	0.088
Total		9300101	729507	100.000

邮箱 : marketing@motifbiotech.com

网址 : <https://www.motifbiotech.com>

地址 : 苏州市高新区浒墅关镇浒青路 111 号 B04