


MS Spectrum



Sample Information

Dissolution method :15%ACN+85%H₂O
Date Acquired :2023/08/22 15:17:26
Injection Volume :0.2ul
Block Temp :200

Interface :ESI
Nebulizing Gas Flow :1.50L/min
CDL Temp :250C
CDL Volt :0v

Prerod Bias :+4.5kv
Detector :-.0.2kv
T.Flow :0.2ml/min
B.conc :50%H₂O/50%MEOH

Name :蓝铜胜肽
Sequence :GHK-Cu (Acetate Salt)
Theoretical :340.37
Observed :340.30