

CERTIFICATE OF ANALYSIS

Product Name	三肽-41
Sequence	FWY (Acetate Salt)
Dissolution condition	5%HCOOH+15%ACN+80%H2O
Length	3AA
Modification	acetate salt
Molecular Weight (MW)	514.58
Storage	-20°C

Test Items	Specifications	Results
MW by MS	514.40	Conforms
Purity by HPLC	>95%	99.796%
Peptide content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	1g	20.0mg*50管

Certified by:

Quality Assurance Department

**Note: this product is intended for research use only; not for diagnostic or human use.**

邮箱: marketing@motifbiotech.com

网址: https://www.motifbiotech.com

地址: 苏州市高新区浒墅关镇浒青路 111 号 B04