


Sample Information

Name	:三肽-41		
Sequence	:FWY (Acetate Salt)		
Pump A	:0.1%Trifluoroacetic in 100% water		
Pump B	:0.1%Trifluoroacetic in 100% acetonitrile		
Total Flow	:1ml/min		
Wavelength	:214nm		
Analytical column type	:Kromasil 100-5-C18 (4.6*250mm*5 µm)		
Dissolution method	:5%HCOOH+15%ACN+80%H2O		
Inj. Volume	: 3ul		
Time	Module	Action	Value
0.01	Pumps	B.Conc	25
30.00	Pumps	B.Conc	55
33.00	Pumps	B.Conc	100
38.00	Pumps	B.Conc	100
40.00	Pumps	B.Conc	25
50.00	Controller	Stop	

Chromatogram



Peak Table

Detector A 214nm

Peak#	Ret. Time	Area	Height	Area%
1	9.638	11509862	1326865	99.796
2	9.904	16133	5396	0.140
3	11.094	7392	1049	0.064
Total		11533387	1333310	100.000

邮箱：marketing@motifbiotech.com

网址：https://www.motifbiotech.com

地址：苏州市高新区浒墅关镇浒青路 111 号 B04