

CERTIFICATE OF ANALYSIS

Product Name	棕榈酰四肽-7
Sequence	Pal-GQPR
Dissolution condition	5%HCOOH+15%ACN+80%H2O
Length	4AA
Modification	HCl Salt
Molecular Weight (MW)	649.92
Storage	-20°C

Test Items	Specifications	Results
MW by MS	649.92	Conforms
Purity by HPLC	>98%	99.10%
Peptide content	>80%	98.18%
Moisture content	< 7.0%	1.18%
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	1g	1.0g

Certified by:
Quality Assurance Department



Note: this product is intended for research use only; not for diagnostic or human use.

邮箱: marketing@motifbiotech.com 网址: <https://www.motifbiotech.com>

地址: 苏州市高新区浒墅关镇浒青路 111 号 B04